



# Welcome to SightSync

Let's Sync In →

Login

Create Account



# Let's Get You Set Up

email

password

Forgot Your Password?

Continue

# SightSync is ready to guide you

Follow the on-screen steps or use voice guidance  
to complete your setup.



Start





Make sure Bluetooth is turned on

---

Select your device

Uzair's Fold Z 5

Malav's Iphone

Pal's Iphone

Shlok's Vivo V60



Voice guidance enabled



Malav's Iphone

Confirm the pairing code.


6 6 8 0 0 1

Continue



Voice guidance enabled





---

## Get to Know Your Button

- single press to read text
- ▬ press and hold to describe your surroundings

---

Setup Complete



Voice guidance enabled



Device

# Malav's SightSync

 Connection Status - Connected

Battery - 85% 

Last Synced - 5 mins ago 





# Controls




Single Press

Describe Scene 



Double Press

Read Text 



Long Press

Off 







# User Guide

Quick Start Guide



Replay Tutorial



How to operate the glasses





# Settings


Audio



Device Management



Account





# Audio

Headphone Volume



Speech Rate



0.5x 1x 1.5x

Voice Selection





# Device Management

Check for Firmware Update



V1.023 - Software Upto Date


Wi-Fi Network



Device Name



Unpair This Device





# Account

Logout



Delete Account

This action cannot be undone

